

# MAGLUMI 800

Chemiluminescence Immunoassay (CLIA) System



# OUTSTANDING TECHNOLOGY POWER OF MAGLUMI

## Chemiluminescence Immunoassay (CLIA) System

CLIA uses two important technologies, one is labeling technology which determines reaction mode; and the other is separation technology which determines the sensitivity, accuracy and precision of the reagents.

### Labeling Technology

Two types of labeling technologies are commonly used. One is enzyme label and the other is non-enzyme small molecule label. Enzyme label reagents are not so stable and are easily affected by the change of storage conditions. Maglumi system applies ABEI labels. ABEI is a non-enzyme small molecule with special molecular formula to enhance stability in acid and alkaline solution. As fast chemiluminescence, ABEI chemical reaction with sodium hydroxide (NaOH) and Hydroperoxide ( $H_2O_2$ ) finishes the process in 3 seconds.

### Separation Technology

Maglumi uses Magnetic Microbeads. As separation technology, it has been widely used in the field of CLIA. Compared with traditional separation technology, it has the following advantages:


- Shortening the reaction time by enlarging the reaction area of antigens and antibodies.
- Enhancing the sensitivity by better and faster capture of antigens and antibodies.
- Reducing inter or intra-assay discrepancies significantly by mixing the reagents thoroughly in a liquid separation platform.
- Enhancing the accuracy by absorbing antigens and antibodies through chemical reaction.



# ADVANTAGE OF MAGLUMI 800



- For small and middle size labs or hospitals
- Simple, Smart, Humanization designed model
- More than enough test menu to meet the different needs and it is enlarging according to your requirement
- Free Quality Control and calibrator helps to reduce the Cost Per Test
- Quality guarantee by the Third party quality control and External Quality Assessment (EQA)
- Key technology ABEI and Nano magnetic microbead enhances stability and sensitivity of MAGLUMI reagents



# MAGLUMI 800

- Integrated design and humanized appearance
- Simple and smart operation
- Compact size and weight



## Pipettor

- Aspirate reagent/sample with high speed
- Titanium needle for crushproof
- Clot detection & liquid level detection
- Coated with TEFLON to prevent carry-over
- Inner & external washing



## Chamber

- High-sensitivity and low-noise PMT
- Cuvettes detection and anti-overflow



## Washer

- Strong magnetic field helps better separation
- Three pairs of washing needles with anti-overflow function



## Automated cuvette loader

- Up to 40 modules of 240 tests
- Continuous loading



## Sample area

- Up to 40 sample on board with barcode
- Continuous loading & cooling function




## Incubator

- Incubate 13 slots (78 tests at the same time)
- The temperature is 36.8 ±0.5



## Reagent area

- Up to 9 reagents on board with RFID
- RFID label storing reagent information
- Easy and fast for adding and updating
- RFID tag, high reliability



# OPERATION SOFTWARE



## Comprehensive software

- User-friendly Interface
- Real-time status monitoring for each test
- Monitoring reagent and consumable status
- Intelligent alarm function



## Maintenance Guide

- Daily, weekly and monthly maintenance guide for analytical unit
- Ensure performance reliability and reduce unnecessary service calls



## Quality control

- Westgard rules and Levey-Jennings chart for both internal and external quality control



## Test Summary Function

- Test summary including system test, calibration, QC, statistics of samples, valid tests and rerun tests
- Search and review test information conveniently



## Pre or post defined dilution Function

- Pre or post defined dilution flexibly if high concentration is met



## Connect to LIS ( bidirection )

- Bidirection, Lis can connect by online and ASTM

# TEST MENU

## Tumor Markers

Ferritin  
AFP  
CEA  
Total PSA  
f-PSA  
CA 125  
CA 15-3  
CA 19-9  
PAP  
CA 50  
CYFRA 21-1  
CA 242  
CA 72-4  
NSE  
S-100  
SCCA  
TPA-snibe  
Pepsinogen I  
Pepsinogen II  
Gastrin-17  
H.pylori IgG

## Bone Metabolism

Intact PTH  
Calcitonin  
Osteocalcin  
25-OH Vitamin D

## Glyco Metabolism

C-Peptide  
Insulin  
ICA  
IAA  
Proinsulin  
GAD 65  
IA-2

## Fertility

FSH  
LH  
HCG/ $\beta$ -HCG  
PRL  
Estradiol  
free Estriol  
Progesterone  
Testosterone  
free Testosterone  
DHEA-S  
**\*AMH**  
17-OH Progesterone

## Thyroid

TSH  
T4  
T3  
FT4  
FT3  
TG  
TGA  
TRAb  
TMA  
Anti-TPO  
Rev T3

## Prenatal Screening

cAFP  
free  $\beta$ -HCG  
PAPP-A

## Inflammation Monitoring

CRP  
PCT

## Kidney Function

$\beta$ 2-MG  
Albumin

## Hepatic Fibrosis

HA  
PIIIP N-P  
C IV  
Laminin  
Cholyglycine

## TORCH

Toxo IgG  
Toxo IgM  
Rubella IgG  
Rubella IgM  
CMV IgG  
CMV IgM  
HSV-1/2 IgG  
HSV-2 IgG  
HSV-1/2 IgM

## EBV

EBV EA IgG  
EBV EA IgA  
EBV VCA IgG  
EBV VCA IgM  
EBV VCA IgA  
EBV NA IgG

## Immunoglobulin

IgM  
IgA  
IgE  
IgG

## Anemia

Vitamin B12  
Ferritin  
FA

## Drug Monitoring

Cyclosporine A  
Tacrolimus, FK 506  
Digoxin

## Cardiac

CK-MB  
Troponin I  
Myoglobin  
NT-proBNP  
Aldosterone  
Angiotensin I  
Angiotensin II  
D-Dimer  
**\*Direct Renin**  
LP-PLA2  
hs-cTnI

## Infectious Disease

HBsAg  
Anti-HBs  
HBeAg  
Anti-HBe  
Anti-HBc  
Anti-HCV  
Syphilis  
Chagas  
HTLV I/II  
Anti-HAV  
HAV IgM  
HIV A b  
HIV p24 Ag  
HIV Ab/Ag combi

## Others

GH  
IGF-I  
Cortisol  
ACTH

\*In Red Available Soon

# MAGLUMI 800

## Chemiluminescence Immunoassay (CLIA) System

### Technical Specifications

#### Reagent Feature

- Flash chemiluminescence label-ABEI, with high sensitivity and long stability
- Nano magnetic microbeads separation, fast and efficient

#### Main Features

- Throughput: up to 180 tests/hour
- 24 hours ready to use
- Time to first result: 17 minutes

#### Modes of Operation

- Random, Batch and STAT

#### Sample Loading

- Up to 40 sample tubes
- Continuous loading, STAT available
- Barcode reader recognition or analyzer auto numbered
- LIS connection, auto read sample info
- Refrigerated sample area with independent power supply

#### Reagent Loading

- 9 reagents on board
- Continuous loading
- RFID reading all info of reagent
- Refrigerated reagent area

#### Reagent Integral

- Ready-to-use, no pretreatment required
- Including calibrators
- RFID tag storing all info of reagent
- RFID tag with built-in master curve
- 2-point calibration to adjust master curve
- Calibration stability: Max 4 weeks

#### Other Features

- Clot detection
- Liquid level detection
- Auto and optional ratio for high sample dilution
- Constant 37°C incubation

#### Operation System

- Windows 7
- DELL compatible Dual-core CPU
- Color touch screen

#### Interconnection

- Bi-communication with LIS via ASTM protocol

#### Dimensions Weight

- Input: AC 110-230V, 50/60 Hz, Max 630W
- Dimension and Weight: 102\*72\*56cm 73kg



Shenzhen New Industries Biomedical Engineering Co., Ltd.(Snibe)

21st Floor, Block A, Building 1, Shenzhen Software Industry Base,  
No. 1008, Keyuan Road, Nanshan District, Shenzhen 518000 CHINA

Tel: +86 755 26501514

Fax: +86 755 26654850

Email: sales@snibe.com

Web: www.snibe.com



DISTRIBUTOR: